



Dryz® Blu post rinse
PHOTO: Dr. F. Shull, SC, USA.

Create a **Dry** Working Environment Every Time.

Dryz®

Gingival Hemostatic Retraction Paste



NEW

- New 'blu' color makes it even easier to see the material.
- Stops gingival bleeding and seepage quickly.
- Aluminum Chloride acts fast and won't slough tissue like Ferric Sulfates.
- New creamy consistency adapts well and washes off clean.
- Available in syringes and single dose capsules, which eliminate waste and keep asepsis.



New Dryz Blu features a bright, contrasting aqua blue color making the material more visible against tooth structure. The slight stiffness of the paste provides better displacement of tissue when the material is pressed subgingivally (with either a compression cap or cotton roll).

The material is effective within 2-3 minutes and rinses easily. Because of the aqua blue color, any residual is quickly spotted and removed simply with a water/air spray.

Ordering Information

Dryz Blu	
S190	● Dryz® Blu Syringes (7) 0.5ml syringes and (15) applicator tips.
S191	● Dryz® Blu Syringes Value Pack (25) 0.5ml syringes and (55) applicator tips.
S181	Dryz® Luer Lock Applicator Tips (40) applicator tips.